VOORSCHRIFTEN:

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**EXPERIMENTELE
IMMUNO-BIOLOGIE**

1. = Preparation of a lymphocyte suspension for lymphocyte cytotoxicity, or frozen storage
2. = Lymphocyten bereiding rat
3. = Mitomycin C treatment of lymphocytes for one-way stimulation in lymphocyte cultures
4. = Oog plexus punctie rat 2 pag.
5. = Selectieve erythrocyten agglutinatie
6. = Preservation of lymphocytes at cryogenic temperatures 2 pag.
7. = Instant tissue culture powder medium
8. = Invriezen tumor celsuspensie
9. = Trypsinatie cellen in suspensie (Lymphocyten of fibroblasten)
10. = Micro cytotoxie ratte lymphocyten
11. = Thymocyten bereiding ratten
12. = Mixed haemadsorptie 2 pag.
13. = Treatment of the cells for DNA-synthesis determination
14. = Lymphocyte culture for determination of ALS activity
15. = Preparation of the lymphocyte suspension for lymphocyte cultures

PREPARATION OF A LYMPHOCYTE SUSPENSION FOR LYMPHOCYTE CYTOTOXICITY,
OR FROZEN STORAGE

- =====
- 1) Venous heparinized blood - heparin without preservative
 "Thromboliquine" (Organon, Oss, Holland)
 8 ml 1:20 in saline
 per 500 ml of blood

 - 2) Removal of the granulocytes by filtration through nylon wool
 either "Leucopak" (Fenwall Laboratory
 Morton Grove, Ill. USA)
 or "Filtralon" (ARA, Amsterdam, Holland)

 - 3) Sedimentation of the red cells at $37^{\circ}\text{C} + 45-60$ min.
 either spontaneously
 or after addition of $\frac{1}{4}$ vol. 5% Dextran
 (MW 180.000 Poviet, Amsterdam)

 - 4) Centrifugation of the lymphocyte rich plasma 10 min 400 G

 - 5) Two washings in Earle's BSS with 10% human AB serum 2 x 10 min
 260 G

 - 6) Resuspension in the media to be used afterwards

 - 7) Cell counting (for quantitative lymphocyte culture experiments
 electronic counting is required)

 - 8) Resuspension to the finally required concentration

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MITOMYCIN C TREATMENT OF LYMPHOCYTES FOR ONE-WAY STIMULATION
IN LYMPHOCYTE CULTURES

=====

The lymphocytes are resuspended in Earle's B.S.S. with mitomycin C (Christiaens, SA, Brussels). Final cell concentration 10×10^6 /ml

Final mitomycin C concentration

12.5 μ g/ml.

Incubation for 30 min. at 37°C.

The cells are spun down, washed once in Earle's B.S.S. (260 x G 10 min), and resuspended in tissue culture medium to a final concentration of 750,000/ml.

In our experiments we culture $1,5 \times 10^6$ mitomycin treated cells together with $1,5 \times 10^6$ non-treated cells from the other donor per tube.

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PRESERVATION OF LYMPHOCYTES AT CRYOGENIC TEMPERATURES
=====

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2. Preparation of the lymphocyte suspension

Venous heparinized blood is used (heparin without preservative e.g. "Thromboliquine" Organon, Oss, Holland) 8 ml 1:20 in saline solution per 500 ml blood. The removal of the granulocytes is accomplished by filtration through "Leukopak" filters (Fenwal Laboratories, Morton Grove, Illinois, USA) at 37°C.

Sedimentation of the red cells is performed at 37°C for 30 - 45 minutes, after addition of 1/4 vol. 5% Dextran (molecular weight 180,000 N.V. Poviet, Amsterdam); then centrifugation of the lymphocyte-rich plasma for 10 minutes at 400 g. followed by two washings in Hanks Balanced Salt Solution for 10 minutes at 400 g. Resuspend the cells in Hanks BSS with 10% normal AB serum and 10% DMSO. The concentration of the cells seems not to be critical. The cell suspension is handled at temperatures of melting ice until the cooling process has started.

The cell suspension is put either in glass sealed 2 ml ampules or in straws holding 1.1 ml. The ampules of these straws are put in racks in the BF3-2 freezer and special care should be taken that enough environmental space is present to guarantee a good heat exchange.

The lymphocytes are cooled at 1°C per minute till -25°C , next 5-7 degrees per minute till -80°C and then the cells are transferred into the liquid nitrogen. In our system, the input of liquid nitrogen is controlled by the Siemens controller and the thermal probe is only used to record that particular cooling rate in the ampule. There is no feed-back from this probe.

The program of our controller is chosen in such a way that a complete or nearly complete compensation of the heat exchange is achieved and that the cooling rates mentioned earlier are maintained. The thawing of the cells is done by immersion in a waterbath at 37°C .

Once the icenucleus has disappeared the suspension is transferred to a tube or a small Erlenmeyer flask and diluted slowly (which is very important) with Hanks BSS with 10% AB serum. After dilution 1:10 with the final concentration of DMO 1%, the cells are spun down and resuspended in either fresh complement or another appropriate vehicle.

INSTANT TISSUE CULTURE POWDER MEDIUM
=====

To prepare IX Liquid Medium

- 1) Measure out 5% less water (triple distilled) than desired total volume of medium. Using a mixing container that is as close to the final volume as possible.
- 2) Add powder medium to 15-20C. (room temperature) water with gentle stirring (do not heat water)
- 3) Rinse out inside of package to remove all traces of powder medium.
- 4) Add 2,2 gms. of NaHCO_3 per liter of medium.
- 5) Dilute to desired volume with water. Stir until dissolved. (do not over-mix)
- 6) Adjust pH of medium to 0.2-0.3 below desired final working pH*; use of 1N NaOH or CO_2 gas is recommended. (add slowly with stirring). After pH has been adjusted, keep container closed until medium is filtered.
- 7) Sterilize immediately by membrane filtration (air pressure system recommended).

* pH units will usually rise 0.1-0.3 upon filtration

If TRIS-HCl buffer is preferred instead of bicarbonate:

Take 15 grams TRIS, add 9.2 cc HCl 36% and dilute with saline to 100 ml. Add 20 ml of this solution to 1000 ml of MEM.
Adjust to final pH with NaOH.

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TRYP SINATIE CELLEN IN SUSPENSIE (LYMPHOCYTEN OF FIBROBLASTEN)
=====

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MICRO CYTOTOXIE RATTE LYMPHOCYTEN

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TREATMENT OF THE CELLS FOR DNA-SYNTHESIS DETERMINATION
=====

24 hours before harvesting - addition of 0.06 μ c thymidine-2-C₁₄
spec. act. 1.2 5 mC/mM per 4 ml culture

- harvesting - one washing in 2% acetic acid
- one washing in saline with 0.5% human AB serum
 - precipitation with 6 ml 10% TCA
 - centrifugation 10' on 3000 rpm
 - removal of the supernatant
 - to 0.6 ml of the sediment 0.3 ml 10% NaOH is added,
whereafter 3 ml cellosolve with 4% formic acid is added
 - this is transferred to a counting vial and the tube is
rinsed with another 3 ml cellosolve formic acid mixture
 - finally 40 ml of scintillation counting fluid is added to
the vial whereafter counting can be performed.

LYMPHOCYTE CULTURE FOR DETERMINATION OF ALS ACTIVITY
=====

Incubation of 0.4 ml cell suspension (containing 3×10^6 lymphocytes) with 0.4 ml of the ALS dilution for 2 hours at 37°C.

Either the cells are washed in excess Earle's BSS whereafter 4 ml of culture medium is added or the cell-serum mixture is diluted without washing with 3.5 ml of culture medium.

PREPARATION OF THE LYMPHOCYTE SUSPENSION FOR LYMPHOCYTE CULTURES
=====

Venous heparinized blood - Heparin without preservative, e.g.
"Thromboliquine" (Organon, OSS, Holland)
16 ml 1:20 in saline solution per 500 ml blood

Removal of the granulocytes by filtration through "Leucopak" filter
(Fenwall Lab., Morton Grove, Ill. USA) at 37°C

Sedimentation of the red cells at 37°C + 45 - 60 min.
either spontaneously
or after addition of 1/4 vol. 5% Dextran
(M.W. 180.000, N.V. Poviet, Amsterdam)

Centrifugation of lymphocyte-rich plasma 10 min. 400 g

Two washings in Hanks B.S.S. 2 x 10 min. 260 g

Resuspension: Minimum Essential Medium (Eagle) for suspension cultures
without L-glutamine (Grand Isl. Biol. Co. N.Y., USA)
+ 20% heat inactivated human AB
+ L-glutamine 0.1 mM/ml
+ penicillin 100 U/ml
+ streptomycin 100 ug/ml

Careful cell counting

Further resuspension in the abovementioned mixture to a final
concentration of 750.000 lymphocytes/ml