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# INSTRUCTIONS FOR YSI 4004 CLARK OXYGEN PROBE

The YSI 4004 Clark Oxygen Probe is a complete polarographic system, comprising a platinum cathode, silver anode, and KCl solution held captive around the electrodes by a Teflon® membrane.

## SPECIFICATIONS:

Electrode Current: At 30°C and 760 mm.  
 Nitrogen (Zero Current): .05 microampere maximum  
 Air: About 2 microamperes  
 Oxygen: About 10 microamperes

## Membrane Temperature

Coefficient: About 4%/°C  
 Polarizing Voltage: 0.8 Volt  
 Response Time: 90% in 10 seconds  
 Steady state in 30 seconds

Oxygen Depletion Rate:  $8 \times 10^{-11}$  grammes oxygen/second/  
 1 microampere

## HOW IT WORKS:

When a polarizing voltage is applied across the probe, all oxygen in the probe is consumed (reduced) at the cathode and current flows in direct stoichiometric relation to the rate of oxygen consumption (reduction). Oxygen then diffuses through the membrane at a mass rate proportional to the oxygen pressure outside the probe, since oxygen pressure in the probe is near zero. When steady state conditions are reached in about 30 seconds, current flows through the probe at a rate in proportion to the external oxygen pressure.

## PROBE FEATURES:

### Isolation

Since the probe is a complete system in itself, it is relatively unaffected by, and does not offset, its external environment. For this reason, the electrode may also be used to measure oxygen in nonconducting liquids or gases. Further, the electrodes are bathed in a known medium and protected from contamination by the membrane. Thus, the probe will measure oxygen in solutions contaminated by ionic reducing agents and reducing (ion consuming) organic matter. The probe is subject

to interference only from low molecular weight reducing gases such as the halogens and hydrogen sulfide.

## Voltage Plateau

The current-oxygen pressure relationship is essentially independent of polarizing voltage within a certain voltage range. Specifically, the output signal shall change less than 3% when the polarizing voltage is lowered from .8 to .65 volts.


## PREPARATION FOR OPERATION:

1. Remove "O" ring and membrane, inspect probe for cleanliness, rinse with DISTILLED water and shake dry.
2. Dilute saturated KCl solution with an equal part of DISTILLED water. (KCl solution available from YSI contains 1-2 drops of Kodak Photo Flo per ounce to insure wetting).
3. Grasp probe in left hand (A), select a membrane and secure one end against side of probe with thumb. Handle the membrane only by the edges. Stretch membrane up, over, and down (B and C). Secure membrane under left forefinger (D). Roll "O" ring over top of probe and into groove to hold membrane in place (E). Sluggish response may be from a slack or wrinkled membrane or foreign matter between membrane and cathode. Wrinkles can be removed by gently pulling the edges of the membrane. Trim off excess membrane near "O" ring (F).
4. Slide the rubber tube down the barrel to expose filler hole. Fill barrel about 3/4 full of half saturated KCl. An eyedropper or syringe is a handy tool for this operation.
5. Inspect the probe visually for trapped air near the sensing end. To remove trapped air, hold probe in one hand, sensing end tilted down, and gently tap barrel near the sensing end. Tapping all around the barrel will help dislodge all trapped air.
6. For simplified membrane mounting procedures we recommend use of the YSI 5339 Probe Service Kit.

## OPERATION:

### Oxygen Depletion

Reading error in liquids due to oxygen depletion in the area

of the membrane will occur unless adequate stirring is provided. Depletion rate is equivalent to  $8 \times 10^{-11}$  grammes oxygen/sec. per sensor current of 1 microampere.

**Membrane Coefficient**

The approximate  $4\%/^{\circ}\text{C}$  membrane coefficient necessitates good temperature control, and temperature equilibration time must be considered when making changes in the setup.

**Cleaning Procedure**


Clean tarnish, etc., off with cotton-tipped swab dipped in 1:1  $\text{NH}_4\text{OH}$  - water solution and flushing the ammonia solution through the probe. Continual noise of several % full scale may occur with time even though the probe and membrane may be in perfect condition. The silver anode may not be making good reference contact with the KCl solution, and re-cleaning with the ammonia solution is recommended. Rinse well with distilled water after each ammonia wash.

**Storage**

Store the probe for short periods submersed in distilled water. For extended periods, clean well and cover with a membrane (no KCl).

**Repair Facilities**

If you are experiencing difficulty with a YSI product, it may be returned to the YSI Customer Service Department for repair, even if the guarantee has expired. YSI maintains complete facilities for prompt servicing of all YSI products.

**REPLACEMENT PARTS AND ACCESSORIES:**

- YSI 5339 Probe Service Kit - Standard Membrane Material Membrane Mounting "O" Rings. "O" Ring Applicator, KCl with Photo Flo, & KCl Dropper Bottle
- YSI 5352 Membrane Kit - 25 (.001" Thk) Standard Membrane
- YSI 5937 High Sensitivity Membrane Kit - 25 (.0005" Thk) Membranes

**GUARANTEE:**

YSI 4004 Clark Oxygen Probes carry a one-year unconditional guarantee on all workmanship and components. Damage through accident, misuse, or tampering will be repaired at a nominal charge, if possible.

**CIRCUIT INFORMATION:**

The following information is offered as a guide for using the YSI 4004 Clark Oxygen Probe. It is not the only way the probe can be used.

**Voltage Supply**

.8 VDC - readily obtained from flashlight cell or mercury battery. If a line operated supply is used the DC must be well filtered and regulated.

**Sensor Load Resistor**


Restrict the voltage across the load resistor to less than 0.10 volt (10-25 mV preferred). The YSI 4004 Sensor will develop about 2.0 microamperes in air at  $30^{\circ}\text{C}$  and 760 mm (10.0 microamperes in oxygen at  $30^{\circ}\text{C}$  and 760 mm).

**AMPLIFIER OR RECORDER INPUT REQUIREMENT:**

The input impedance of the amplifier should be large compared to the sensor load resistor; if not, the input resistance shunts the sensor load resistor. In many cases this is permissible if the effect is recognized and proper circuit modification made. For instance if the amplifier input impedance is purely resistive and has a value of 100,000 ohms the sensor load resistor may be increased so that the resultant paralleled resistors are the desired value.

NOTE: Wire wound resistors and potentiometers are best - 10-turn potentiometers give better resolution - A switch may be installed to permit easy change from air to oxygen.

Values of load resistor for air and oxygen are for use at temperatures between 30 and  $40^{\circ}\text{C}$ . The signal voltage will be 10 millivolts. If the amplifier requires a larger input voltage for full scale these resistors may be increased proportionately.

**Fig. 1**

4004 Oxygen Electrode available only directly from

**YELLOW SPRINGS INSTRUMENT CO.**  
**YELLOW SPRINGS, OHIO 45387**

Fig. 2 Oxygen measurement chamber